
KPL SignalOCK™ ChemiWestern Kits

(Film and Imager Analysis)

Products	Catalog No.
KPL SignalOCK™ HRP ChemiWestern Kit (Film)	5410-0010 (54-53-00)
KPL SignalOCK™ HRP ChemiWestern Kit (Imager)	5410-0011 (54-54-00)
KPL SignalOCK™ AP ChemiWestern Kit (Film/Imager)	5410-0012 (54-56-00)



TABLE OF CONTENTS

Section	Page
INTRODUCTION.....	4
PRINCIPLE OF KPL SignalLOCK™ ChemiWestern KITS.....	4
MATERIALS AND EQUIPMENT	4
PRODUCT SAFETY AND HANDLING.....	5
REQUIRED SUPPLIES AND EQUIPMENT FOR WESTERN BLOTTING NOT INCLUDED.....	5
WARNINGS AND PRECAUTIONS.....	5
WESTERN BLOT DETECTION AT A GLANCE	5
NOTES FOR OPTIMAL KIT PERFORMANCE	6
POLYACRYLAMIDE GEL ELECTROPHORESIS AND WESTERN BLOTTING.....	6
DETECTION OF WESTERN BLOTS.....	6
STRIPPING AND REPROBING MEMBRANES	8
TROUBLESHOOTING GUIDE FOR CHEMILUMINESCENT.....	8
DETECTION OF WESTERN BLOTS.....	8
PROBLEM 1: NO SIGNAL.....	8
PROBLEM 2: WEAK SIGNAL.....	8
PROBLEM 3: EXCESSIVE SIGNAL, NONSPECIFIC BANDS OR GENERAL BACKGROUND....	9
PROBLEM 4: POORLY DEFINED OR “FUZZY” BANDS OR DOTS	9
RELATED PRODUCTS	9

Introduction

Western blotting is a highly specific method for analysis of proteins. Use of an alkaline phosphatase (AP) or horseradish peroxidase (HRP)-linked affinity purified antibody together with a highly sensitive chemiluminescent substrate provides an excellent method for detection and characterization of samples bound to membranes through Western or dot blotting. Proteins are immobilized on membranes by way of either electrophoresis and subsequent transfer or dot blotting. For dot blots, proteins are spotted directly to the membrane. For a Western blot, proteins are separated by SDS-polyacrylamide gel electrophoresis (or a comparable method) and transferred to the membrane through either electrophoretic or passive transfer. Following attachment to membrane, a primary antibody is used to selectively bind the protein of interest. Alternatively, a known protein is bound to the membrane for screening of specific monoclonal antibodies or serum samples. An enzyme-labeled secondary antibody directed against the species of origin of the primary antibody or antiserum is then applied. This antibody is coupled to an enzyme (either AP or HRP) which work to destabilize the acridane present in each substrate. This reaction allows for the emittance of light. The resulting signal can be measured using a camera or X-ray film.

Principle of KPL SignalLOCK™ ChemiWestern Kits

SeraCare offers a line of elite chemiluminescent substrates, KPL LumiGLO Reserve™, KPL LumiGLO Ultra™ and KPL PhosphaGLO Reserve™ that have been paired with a superior blocker, KPL SignalLOCK to achieve the ultimate in Western blot sensitivity. Femtogram detection is readily achievable with these kits. They offer improvements in the way of signal intensity – greater than 20-fold more sensitive than LumiGLO® and other competitive substrates. These products are specifically designed for the detection of proteins that are (1) difficult to detect because they are in such low quantities; or (2) from samples that are precious and, therefore, desired to be conserved. In addition, these kits provide the benefit of strong signal with the use of reduced amounts of antibodies and antibody conjugates.

Kits have been carefully optimized to allow for ease of selection for either imager or film detection. KPL SignalLOCK HRP ChemiWestern Kit (Film) contains KPL LumiGLO Reserve, a known high-performance substrate, which will accomplish femtogram detection on film. KPL SignalLOCK HRP ChemiWestern Kit (Imager) utilizes KPL LumiGLO Ultra, a novel substrate which allows for rapid, sensitive detection on imagers. KPL SignalLOCK AP ChemiWestern Kit (Film/Imagers) allows for versatile detection on both film and imagers with the use of KPL PhosphaGLO Reserve. All three kits utilize KPL SignalLOCK Blocking Solution, a non-protein blocker that eliminates cross-reactivity issues as well as enhancing signal visibility. For added convenience, KPL Wash Solution is also included to further minimize background issues.

These combinations deliver rapid and accurate identification of proteins for your more challenging detection work where one is forced to manage signal:noise contrast when pushing the detection limits of Western blot to clearly detect low abundance, short half-life, poorly expressed or otherwise small samples of protein.

MATERIALS AND EQUIPMENT

KPL SignalLOCK™ HRP ChemiWestern Kit (Film)

KPL LumiGLO Reserve™ Chemiluminescent Substrate A
KPL LumiGLO Reserve Chemiluminescent Substrate B
KPL 5X SignalLOCK™ Blocking Solution
KPL 20X Wash Solution

KPL SignalLOCK™ HRP ChemiWestern Kit (Imager)

KPL LumiGLO Ultra™ Chemiluminescent Substrate A
KPL LumiGLO Ultra Chemiluminescent Substrate B
KPL 5X SignalLOCK™ Blocking Solution
KPL 20X Wash Solution

KPL SignalLOCK™ AP ChemiWestern Kit (Film/Imager)

KPL PhosphaGLO Reserve Chemiluminescent Substrate
KPL 5X SignalLOCK™ Blocking Solution
KPL 20X Wash Solution Concentrate

Sufficient reagents are provided to test approximately 1000 cm² of membrane (ten 10 cm x 10 cm mini-blots) when using recommended minimal volumes. Reagents are stable for a minimum of one year when stored at 2–8°C. KPL SignalLOCK Blocking Solution is stable for a minimum of one year when stored at 2–25°C.

PRODUCT SAFETY AND HANDLING

See SDS (Material Safety Data Sheet) for this product.

REQUIRED SUPPLIES AND EQUIPMENT FOR WESTERN BLOTTING NOT INCLUDED

- Primary and secondary antibodies (purified)
- Nitrocellulose or PVDF Membrane (Nylon is not recommended)
- Platform shaker or rocker
- Blotting or Whatman paper
- Electrophoresis equipment
- Polyacrylamide gels
- Protein standards
- X-ray film (double emulsion) or CCD Imager
- Developing chemicals/equipment
- Incubation trays or tubes (Note: Gel box containers are not suitable for incubating membranes).

WARNINGS AND PRECAUTIONS

- Read ALL instructions thoroughly before using the kit.
- Always wear protective gloves and a lab coat for personal protection, as well as protection of the membrane and immunoassay reagents from contaminants such as skin oils or proteins.
- Caution should be used when preparing or handling polyacrylamide gels; monomeric acrylamide is a neurotoxin.
- For proper analysis of results, always include positive and negative controls, blanks and/or protein standards as appropriate.
- Prior to application of the kit reagents, the protein of interest must be immobilized onto the test membrane. Nitrocellulose and polyvinylidene difluoride (PVDF) membranes have both been determined to be suitable for use with these reagents.

WESTERN BLOT DETECTION AT A GLANCE

Total Time: 4 hours

Polyacrylamide Gel Electrophoresis



Immobilize Protein on Membrane



Block Membrane

1 hour or overnight



Incubate Primary Antibody

30 minutes - 1 hour



Wash Membrane

3 x 5 minutes per wash



Incubate Conjugate

30 minutes - 1 hour



Wash Membrane

1X Wash Solution, 4 x 5 minutes per wash

Reagent Quality Water, 1 x 5 minutes



Incubate with Substrate

1 minute



Expose to Film or Imager

10 seconds - 10 minutes

NOTES FOR OPTIMAL KIT PERFORMANCE

- To achieve optimal results, please blot the membrane completely dry with filter paper before image analysis.
- Different kits require different ratios of substrate. Special instructions for mixing of substrates are highlighted as **IMPORTANT: Substrate Instructions**.
- It is imperative to optimize primary and conjugate dilutions for your Western blot assay. Slight differences in activity can result in overwhelming background when employing powerful chemiluminescent substrates for detection. Primary and secondary antibodies should be purified to ensure optimal results.

Polyacrylamide Gel Electrophoresis and Western Blotting

STEPS	CRITICAL POINTS
1. Prepare samples by diluting to desired concentration with sample diluent. If using a reducing agent (ex. – β - mercaptoethanol), incubate samples at 100°C for 3 minutes prior to electrophoresis.	<i>If this kit is being used for the first time, reduced volumes of target may be loaded to achieve comparable detection to routine HRP-chemiluminescent substrates.</i>
2. Electrophorese samples and transfer by standard methods.	<i>Both nitrocellulose and PVDF membranes may be used. Nitrocellulose is preferred.</i>
3. Mark the orientation of the protein samples on the membrane prior to detection.	<i>The membrane may be cut into strips at this time if desired. Alternatively, it may be more convenient to cut strips after the entire membrane has been blocked.</i>

Detection of Western Blots

There are many protocols available for the detection of Western blots. For optimal signal to noise and sensitivity, the following protocol and reagents are recommended.

STEPS	CRITICAL POINTS
1. Block the membrane by immersing in 1X KPL SignalOCK using a minimum of 0.2 mL/cm ² of membrane. Block at room temperature for 1 hour with gentle rocking or shaking, or stationary at 2–8°C overnight. If blocking overnight, container should be covered.	<i>Dilute KPL SignalOCK 1:5 in reagent quality water. Make sure to use a container of proper size that allows the block solution to freely float over the membrane.</i>
2. Remove blocking solution. Incubate membrane with purified primary antibody diluted 1:1,000 – 1:10,000 from a 1 mg/mL stock for 30 minutes - 1 hour. This antibody should be diluted in fresh blocking solution.	<i>It is recommended that serial dilutions through a dot blot are performed to determine the optimal working dilution of the primary antibody.</i>
3. Remove primary antibody in block solution. Wash the membrane in a generous amount of 1X KPL Wash Solution (at least 25 mL for a 100 cm ² membrane). Wash membrane 3 times for 5 minutes per wash with agitation.	<i>Dilute 20X KPL Wash Solution 1:20 in reagent quality water. This solution will provide optimal signal to noise.</i>
4. Dilute appropriate enzyme-labeled secondary antibody 1/25,000 – 1/250,000 of a 1.0 mg/mL stock in freshly prepared 1X KPL SignalOCK using a minimum of 0.2 mL/cm ² of membrane.	<i>The optimal dilution may vary for different lots of conjugate. It is imperative that you titrate the conjugate to determine the optimal working dilution. These dilutions hold true for either AP or HRP.</i>
5. Incubate blot with diluted conjugate for 30 minutes - 1 hour at room temperature with shaking.	
6. Remove the substrate from refrigeration. During the conjugate incubation step, prepare the substrate.	<i>Allow substrate to warm to room temperature prior to use. Minimize exposure to light.</i>

STRIPPING AND REPROBING MEMBRANES

After performing protein transfer, detection with substrate and film exposure, membranes may be stripped and reprobed with new primary and secondary antibodies.

1. Strip antibodies by incubating blot for 30 - 90 minutes at 70°C in erasure buffer: 2% SDS (w/v), 62.5 mM Tris-HCl (pH 6.8 at 20°C), 100 mM β-mercaptoethanol.
2. Wash 2 times, for 10 minutes each, in TBS: 10 mM Tris-HCl (pH 7.4 at 20°C), 150 mM NaCl.
3. Block for 2.5 hours in 1X KPL SignalLOCK Block (or equivalent block).
4. Repeat detection procedure.

TROUBLESHOOTING GUIDE FOR CHEMILUMINESCENT DETECTION OF WESTERN BLOTS

PROBLEM 1: NO SIGNAL

Possible Cause	Corrective Measure
<ul style="list-style-type: none">• Inactive enzyme conjugate.	
<ul style="list-style-type: none">• No binding of conjugate to the primary antibody.	

PROBLEM 3: EXCESSIVE SIGNAL, NONSPECIFIC BANDS OR GENERAL BACKGROUND

Possible Cause

Corrective Measure

- | | |
|---|--|
| <ul style="list-style-type: none">• Excessive antibody used for detection | <p>Optimize conjugate concentration. Reduce antibody concentrations; optimal conjugate dilution should be 1/10,000 – 1/100,000 of a 0.1 mg/mL stock.</p> <p>OR</p> <p>Decrease the amount of primary antibody.</p> |
| <ul style="list-style-type: none">• Overexposure of blot | <p>1</p> <p>2</p> |
| <ul style="list-style-type: none">• Excessive protein loaded on the gel | |
| <ul style="list-style-type: none">• Insufficient blocking or washing, causing non-specific reaction | |
| <ul style="list-style-type: none">• Endogenous peroxidase in the sample (HRP only) | |
| <ul style="list-style-type: none">• Extra bands not due to blocker. | |

KPL LumiGLO is a registered trademark of SeraCare Life Sciences, Inc.

KPL SignalOCK, KPL ReserveAP, KPL PhosphaGLO Reserve, KPL LumiGLO Ultra and KPL LumiGLO Reserve are trademarks of SeraCare Life Sciences, Inc.

TWEEN is a trademark of ICI Americas Inc.

Note: The recommendations of this bulletin are provided solely for the benefit of users who need practical guidance on immunoassay procedures. Due to the fact that experimental conditions for the use of the suggested products are beyond the control of SeraCare Life Sciences, Inc., it is impossible for SeraCare to implicitly guarantee the performance of the mentioned products for any and all assay procedures. Users who need additional information are encouraged to call Technical Services at 800/638-3167 or 301/591-8200 for assistance.

