

# KPL BacTrace® Anti-*E. coli* O111 Magnetic Beads

Catalog No.  
5350-0039 (082-01-95-91)

Size  
1 mL

## DESCRIPTION

KPL BacTrace Anti-*Escherichia coli* O111 Magnetic Beads are super-paramagnetic polystyrene beads coated with SeraCare's KPL BacTrace Anti-*E. coli* O111 antibody. They are intended for use in the isolation and separation of *E. coli* O111 from a variety of samples (food, animal feed, environmental samples, etc.). Immunomagnetic separation offers a rapid means of separating *E. coli* O111 from complex mixtures prior to immunodetection assays, PCR or other culture techniques.

## FORM/STORAGE

Suspension. Store at 2-8°C. DO NOT FREEZE! Stable for a minimum of 1 year from date of receipt when stored at 2-8°C.

## STABILIZER AND PRESERVATIVE

Bovine Serum Albumin (BSA) is added as a stabilizer. 0.02% sodium azide is added as a preservative. Non-sterile.

## BEAD CONCENTRATION

Beads are provided at a concentration of  $> 1 \times 10^9$  beads/mL. Beads are approximately 2  $\mu\text{m}$  in size. One vial provides enough material to perform 50 extractions when using 1 mL of enriched culture.

## MATERIALS REQUIRED BUT NOT PROVIDED

### Wash Buffer:

- 0.5 g BSA, 50  $\mu\text{L}$  Tween 20 in 100 mL Buffered Peptone Water
- Magnetic Separator
- Sterile 1.5 mL Microfuge Tubes
- 1 mL Pipette and Sterile Tips
- 20  $\mu\text{L}$  – 200  $\mu\text{L}$  Pipette and Sterile Tips
- Lab Rotator

## SUGGESTED PROTOCOL

**PLEASE NOTE: Working with pathogenic bacteria requires that certain safety measures be followed. Please follow all required aseptic techniques, as well as good laboratory practice. Endeavor to avoid aerosol formation, and perform necessary work in a biosafety cabinet. All contaminated materials should be autoclaved or disinfected prior to disposal. Follow all pertinent regulations.**

**Prior to use:** perform enrichment steps according to established protocols.

1. After enrichment, homogenize the sample as recommended (e.g. – stomacher). Allow sample to settle for 2 – 5 minutes.
2. Resuspend magnetic beads by inverting several times or vortexing.
3. Pipette 1 mL of homogenized culture into a sterile microfuge tube, taking care to avoid any debris remaining in the sample.
4. Add 20  $\mu\text{L}$  of magnetic beads to the sample. Incubate while rotating for 15 minutes.
5. Place the tube in a magnetic separator for 3 minutes.
6. Carefully remove the supernatant from the tube and discard.
7. Add 1 mL of wash buffer.
8. Remove from separator and invert several times.
9. Repeat steps 5-8 four additional times, for a total of five washes.
10. The beads binding bacteria are now ready for plating or other protocols.

Total protocol time is approximately 45 minutes.

**KPL BacTrace®**  
**Anti-*E. coli* O111**  
**Magnetic Beads**

<u>Catalog No.</u>	<u>Size</u>
5350-0039 (082-01-95-91)	1 mL

**REFERENCES**

1. Detection and Isolation of non-O157 Shiga Toxin-Producing *Escherichia coli* (STEC) from Meat Products. USDA SOP No. MLG 5B.03. 2012, pgs. 12 – 15.
2. Gehring, et. al. Enzyme-linked Immunomagnetic Electrochemical Detection of *Salmonella Typhimurium*. *J Immunol Methods*. 1996 Sep 9;195(1-2):15-25.